


# AC440 Plug-in for AuditBase System



The AC440 plug-in from Interacoustics appears as a small control panel that runs on top of the AuditBase Audiogram. The clinician controls the audiometer via the panel, and the result of the measurement is entered directly into the AuditBase Audiogram.

## AC440 Plug-in for AuditBase

The AC440 Plug-in for AuditBase System gives condensed access to the new Interacoustics Affinity<sup>2.0</sup>/Equinox<sup>2.0</sup> AC440 module. Supports the daily operation in clinics that use AuditBase.

### Features:

- Full 2 channel clinical audiometry
- Perform Tone, Speech and Weber tests
- Talk forward, Talk back and Monitor functions
- Extended Range (+20 db) button
- Presentation control: Manual/Reverse, Pulsation, Synchronous, Simultaneous/Alternate.
- Choose between 1 or 5 dB steps
- Select HCL, MCL or UCL
- Advanced speech testing:
  - Manual correct / incorrect scoring
  - SRT, UCL, WR1 & WR2
  - Speech from hard drive and easy word presentation
- Full keyboard control using shortcuts



**Interacoustics**<sup>®</sup>

leading diagnostic solutions

# AC440 Plug-in for AuditBase System

## Specifications:

The AC440 Plug-in for AuditBase System works with the following Interacoustics instruments:


Equinox<sup>2.0</sup> PC-based Clinical Audiometer  
Affinity<sup>2.0</sup> PC-based Hearing Aid Analyzer

## Software Requirements:

The AC440<sup>2.0</sup> module (Included in the Equinox or Affinity Suite)  
AuditBase version 4.14.XX and later.



The AC440 plug-in from Interacoustics appears as a small control panel that runs on top of the AuditBase Audiogram module. The clinician controls the audiometer via the panel, and the result of the measurement is entered directly into the AuditBase Audiogram.



## AuditBase System

AuditBase System is developed, distributed and supported by Auditdata A/S. Visit [www.auditdata.com](http://www.auditdata.com) for more information.

Sales and service in your area:

Read more here:  
[www.interacoustics.com](http://www.interacoustics.com)

Interacoustics A/S

Phone: +45 6371 3555 · Fax: +45 6371 3522

E-mail: [info@interacoustics.com](mailto:info@interacoustics.com)

DK-5610 Assens, Denmark

Web: [www.interacoustics.com](http://www.interacoustics.com)



# Interacoustics®

leading diagnostic solutions